

ONE STEP CAN MAKE THE DIFFERENCE

With your help, a step onto a scale can pave the way to a healthier future. To combat chronic malnutrition and address the root causes of stunting, accurate weight measurement is essential. The seca 874+ is your key to delivering fundamental healthcare in even the most remote areas, supporting your mission to enhance living conditions in underserved regions.

Helping you help.



PRECISION FOR LAST-MILE HEALTHCARE



seca 874+

Recognizing the crucial role of health workers and the challenges they face, we introduce the seca 874+, a mobile scale designed to measure the fundamental health parameter of weight with utmost accuracy, ease, and reliability.

Experience the seca 874+ in the Field



Recon, a non-profit aid organization, is dedicated to improving living conditions in underserved rural communities in Colombia. We provided one of our seca 874+ devices to their health workers for a demonstration in Guajira, a department located in northeastern Colombia.



Watch the video to discover more about this project.



"Drawing on decades of experience as a close partner of aid organizations, we have developed the ideal scale designed to empower community health workers in their vital mission to improve public health."

Antonio Cordova
seca Head of Sales International

The #1 Recommended Scale for Community Health Programs



1. Antimicrobial ABS plastic for better hygiene: Incorporating an antimicrobial additive into the non-slip plastic tread surface inhibits germ spread, ensuring impeccable hygiene.

2. Weight optimization: The metal frame is designed to maximize weight efficiency, providing an ideal balance of lightness and strength for ultimate precision and comfort.

3. Glow-in-the-dark LCD display: The long-lasting luminescent LCD allows perfect weight readings even in low light.

4. 25,000 weighings with only 2 batteries: With its energy-efficient design, the scale is powered by just 2 AA batteries, minimizing your battery consumption compared to conventional scales powered by 6 AA batteries.

5. Optional Bluetooth interface: Connect the seca 874+ easily to your smartphone via Bluetooth and immediately send measurements to a documentation app.

Further Benefits:

Resistant to Mud, Dust and Water



With IP53 protection against dust and splash water, the seca 874+ is durable and dependable, even outdoors.

Weigh Toddlers and Babies Comfortably



Easily weigh babies and toddlers. Let the parent step on the scale, press the foot switch, and then re-weigh with the child.

Leveling Feet for Perfect Stance



Level the scale for optimal precision on rough surfaces. The leveling feet are fixed and can't get lost.

Single Foot Switch & Dual-Sided Display



Easy to use: read the results from either side, with the LCD digits adjustable towards or away from the patient.

Select Your seca 874+ and Find the Perfect Model for Your Needs



seca 874+

seca 874+

Bluetooth

Capacity	250 kg or 200 kg	250 kg or 200 kg
Graduation	100 g	100 g
Dimensions (WxHxD)	360 × 61 × 400 mm / 14 x 2.4 × 15.7 in	360 × 61 × 400 mm / 14 x 2.4 × 15.7 in
Platform (WxHxD)	350 × 57.4 × 315 mm / 13.8 × 2.3 × 12.4 in	350 × 57.4 × 315 mm / 13.8 × 2.3 × 12.4 in
Net weight	3.7 kg / 8.2 lbs	3.7 kg / 8.2 lbs
Power supply	Batteries, Power adapter (optional)	Batteries, Power adapter (optional)
Functions	TARE, Auto-HOLD, CLEAR, Mother+child function, Standby mode	TARE, Auto-HOLD, CLEAR, Mother+child function, Standby mode
Interfaces		Bluetooth®
Accessories	+ Switch-mode power adapter seca 400 + Power adapter 401	+ Switch-mode power adapter seca 400 + Power adapter 401
Medical Device Regulation (MDR) certified		



Since the 1980s, we have been developing scales for aid organizations, giving us a precise understanding of their requirements. The seca 874+ represents the pinnacle of this evolution, incorporating decades of experience and technological advancements.

seca 890



seca 881



seca 874



seca 874+



seca.com