

seca 310

Mechanical circular dial
scales for use far from any
hospital



- Fine graduation of either 1 oz / 25 g, 2 oz / 50 g, 25 g or 50 g.
- Tough and reliable.
- Suspension hooks made of galvanized steel, extremely durable and safe.
- One sling seat seca 410 included in delivery

seca 310:

The flexible companion for traveling health providers

Seca knows what really counts when babies and toddlers have to be weighed: the child's safety and precise results. And in baby's first few months it is very important to keep an accurate check on his or her weight. To make that an easy option at any time and in any place, seca has developed an extremely versatile, precise scale. The suspended scale seca 310 with dial and pointer. It can be used at any location, weighs very little itself, and is extremely robust - exactly right for use outside the clinical environment, in particular by visiting midwives and health workers.



Very versatile

The seca 310 is ready for use at any time. It only needs to be hung up by the safety hook made of high-grade galvanized steel. If there is nowhere to hang the scale, the hook doubles as a convenient handle, so the scale can be held with one hand. To weigh the child, he/she is placed in the included sling seat seca 410 fixed to the lower edge of the scale by a second safety hook.



Small units for highly accurate readings

The dial is 7 inch / 180 mm in diameter and the pointer is precisely aligned – a combination that produces especially clear results. Moreover the dial is so large that small units can be marked, guaranteeing very accurate readings of the child's weight. So it is easy to register exactly even the slightest of changes – a major advantage when monitoring the weight of babies. And the seca suspended scale with dial and pointer is the only one on the market to offer this benefit.

Extremely sturdy design

A scale designed to be carried about must withstand a lot of wear and tear. Which is why the seca 310 is housed in an extremely robust metal casing. A guard ring prevents damage to the glass front of the dial. The glass is flat, which offers especially good protection from scratching. To guarantee that the scale is long-lasting, the metal parts are treated to guard against corrosion and the plastic used is UV resistant.



Optional: carrying case seca 411

Infants sit or lie for no-hassle weighing

The adjustable suspension system and the ergonomic shape make light work of weighing babies (lying down) and toddlers (seated). The sling is very flexible, adjusting perfectly to fit the child's size. The cloth has a netting weave, allowing the health worker to keep an eye on the baby lying in the sling throughout the weighing process. When toddlers are weighed sitting down, the holes for the legs hold them safely in the sling.

Easy to carry and clean

Weighing only 3.4 lbs / 1.5 kg (scale + hook) the seca 310 is a lightweight to carry. Scale and cloth are easy to clean, the nylon sling can be washed on the hottest cycle. A sling seat seca 410 is included with the scale.

seca 310

Technical data

- Capacity:
33 lbs / 15 kg (Type 1), 55 lbs / 25 kg (Type 2),
15 kg (Type 3), 25 kg (Type 4)
- Graduation:
1 oz / 25 g (Type 1), 2 oz / 50 g (Type 2),
25 g (Type 3), 50 g (Type 4)
- Dimensions (WxHxD):
8.7 x 8.9 x 2" / 220 x 225 x 52 mm
- Weight: 3.4 lbs / 1.5 kg
- Functions: reset-to-zero function
- Optional: set of 5 pieces sling seat
seca 410, carrying case seca 411